

Mobility Scooter Battery Charging Instructions

1. Before First Use

- Fully charge the batteries **before using the scooter for the first time**.
 - Initial charge can take **8–12 hours** (or overnight).
 - Must be left on until charging has finished, indicated on the charger.
-

2. How to Charge the Scooter

1. Turn the scooter **OFF**.
 2. Locate the **charging port** (usually on the tiller/dashboard or battery pack).
 3. Plug the **charger into the scooter first**.
 4. Then plug the charger into a **standard wall outlet**.
 5. Turn on the power at the wall, only if the charger has no on/off switch.
-

3. Charging Time

- Typical charge time: **6–10 hours** depending on battery level.
 - When the charger light turns **green**, the batteries are fully charged.
-

4. Daily Use Charging

- Charge the scooter **after every use**, even if only used for a short trip.
 - Even after a short trip you need to recharge your scooter.
 - Batteries can take a while to recharge even after a short trip 1-2 hours
 - A long trip the batteries can take 10-14 hours
 - **Try to avoid short charging times, always leave on charge until charger has finished charging, this can lead to shorter life span for the batteries.**
 - **Once charging is complete, turn off the charger and unplug the charger from the wall.**
-

5. Important Battery Care Tips

- Do **not** leave batteries in a discharged state, this can reduce the life of the battery.
 - Do not let batteries discharge below a $\frac{1}{4}$ on the battery gauge as doing this may reduce the life of the battery.
 - Keep the scooter in a **dry, ventilated area** while charging.
 - Use **only the charger supplied** with the scooter.
-

6. Storage Charging

- Charge the batteries **at least once every 1–2 weeks**.
- Store fully charged, never empty.